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Einstellungen alte Firmware:

DAP-1665 // AP	SETUP	ADVANCED	MAINTENANCE	STATUS	HELP
SETUP WIZARD WIRELESS SETUP LAN SETUP	<div style="background-color: #f4a460; padding: 5px;">WIRELESS NETWORK</div> <p>Use this section to configure the wireless settings for your D-Link Access Point. Please note that changes made on this section may also need to be duplicated on your wireless client.</p> <p>Save Settings Don't Save Settings</p>				Helpful Hints.. Wireless Mode : Select a function mode to configure your wireless network. Function modes include Access Point, AP Client, Bridge, Bridge with AP and Repeater. Function wireless modes are designed to support various wireless network topologies and applications. Wireless Network Name : Changing your Wireless Network Name is the first step in securing your wireless network. We recommend that you change it to a familiar name that does not contain any personal information. Hidden Wireless : Enabling Hidden Mode is another way to secure your network. With this option enabled, no wireless clients will be able to see your wireless network when they perform scan to see what's available. In order for your wireless devices to connect to your AP, you will need to manually enter the Wireless Network Name on each device. Security Keys : If you have enabled Wireless Security, make sure you write down WEP Key or Passphrase that you have configured. You will need to enter this information on any wireless device that you connect to your wireless network. Bridge setting : If you want to bridge with the other Bridge AP, please write down the MAC address of the Bridge AP. Besides, you also need to write down the MAC
	WIRELESS MODE : Wireless Mode : Access Point Site Survey				
	2.4GHZ WIRELESS NETWORK SETTINGS : Enable Wireless : <input checked="" type="checkbox"/> Always Add New Wireless Network Name : XRAY_SWISS (Also called the SSID) 802.11 Mode : Mixed 802.11n, 802.11g and 802.11b Wireless Channel : 6 Enable Auto Channel Scan : <input checked="" type="checkbox"/> Channel Width : Auto 20/40MHz Visibility Status : <input checked="" type="radio"/> Visible <input type="radio"/> Invisible				
	2.4GHZ WIRELESS SECURITY SETTING : Security Mode : WPA-Personal WPA Mode : AUTO (WPA or WPA2) Cipher Type : TKIP and AES Pre-Shared Key : ••••••••				
	5GHZ WIRELESS NETWORK SETTINGS : Enable Wireless : <input checked="" type="checkbox"/> Always Add New Wireless Network Name : XRAY_SWISS (Also called the SSID) 802.11 Mode : Mixed 802.11ac, 802.11n and 802.11a Wireless Channel : 36 Enable Auto Channel Scan : <input checked="" type="checkbox"/> Channel Width : Auto 20/40/80 MHz Visibility Status : <input checked="" type="radio"/> Visible <input type="radio"/> Invisible				
	5GHZ WIRELESS SECURITY SETTING : Security Mode : WPA-Personal WPA Mode : AUTO (WPA or WPA2) Cipher Type : TKIP and AES Pre-Shared Key : ••••••••				

Einstellungen alte Firmware:

DAP-1665 HW:B1 FW:2.00

[Home](#) | [Settings](#) | [Management](#)

Existing Network Connected

Click on any item in the diagram for more information.

DAP-1665

IPv4 Network

MAC Address: 80.26.89.B3.CA.10
 IP Address: 192.168.22.3
 Subnet Mask: 255.255.255.0
 Default Gateway: Not Available
 Primary DNS Server: Not Available
 Secondary DNS Server: Not Available

IPv6 Network

Link-Local Address: FE80::8226:89FF:FEB3:CA10
 Router IPv6 Address: Not Available
 Default Gateway: Not Available
 Primary DNS Server: Not Available
 Secondary DNS Server: Not Available

Wi-Fi 2.4GHz

Status: Disabled
 Wi-Fi Name (SSID): Not Available
 Password: Not Available

Wi-Fi 5GHz

Status: Enabled
 Wi-Fi Name (SSID): DIGIPAX
 Password: avmx05405

[Go to settings](#)

[Go to settings](#)

Extender

Use this section to configure the wireless settings for your D-Link Ap. Please note that changes made on this section may also need to be duplicated on your Wireless Client.

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The screenshot displays the D-Link DAP-2660 web management interface. The top navigation bar includes 'Home', 'Maintenance', 'Configuration', 'System', 'Logout', and 'Help'. A left sidebar shows a tree view with 'DAP-2660' expanded to 'Basic Settings', which includes 'Wireless', 'LAN', and 'IPv6'. Below this are 'Advanced Settings' and 'Status'. The main content area is titled 'Administration Settings' and contains several sections, each with a checkbox: 'Limit Administrator', 'System Name Settings', 'Login Settings' (checked), 'Console Settings', 'SNMP Settings', 'Ping Control Setting', 'LED Settings', and 'Central WiFiManager Setting'. The 'Login Settings' section is active, showing 'Login Name' as 'admin', 'New Password' as 'digi4admin', and an empty 'Confirm Password' field. There is an 'Apply New Password' checkbox. A 'Save' button is located at the bottom right of the main content area.

LAN Konfiguration Router:

The screenshot displays the D-Link DAP-2660 web management interface. The top navigation bar includes 'Home', 'Maintenance', 'Configuration', 'System', 'Logout', and 'Help'. A sidebar on the left shows a tree view with 'DAP-2660' expanded, containing 'Basic Settings' (Wireless, LAN, IPv6), 'Advanced Settings', and 'Status'. The main content area is titled 'LAN Settings' and contains the following configuration fields:

Get IP From	Static IP (Manual)
IP Address	192.168.22.3
Subnet Mask	255.255.255.0
Default Gateway	
DNS	

A 'Save' button is located at the bottom right of the configuration area.

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The screenshot shows the D-Link DAP-2660 configuration interface. The main menu includes Home, Maintenance, Configuration, System, Logout, and Help. The left sidebar shows a tree view with DAP-2660, Basic Settings (Wireless, LAN, IPv6), Advanced Settings, and Status. The main content area is titled 'Wireless Settings' and contains the following fields:

- Wireless Band: 5GHz
- Mode: Access Point
- Network Name (SSID): DIGIPAX
- SSID Visibility: Enable
- Auto Channel Selection: Enable
- Channel: 100
- Channel Width: Auto 20/40/80 MHz
- Authentication: WPA-Personal
- PassPhrase Settings:
 - WPA Mode: AUTO (WPA or WPA2)
 - Cipher Type: Auto
 - Group Key Update Interval: 3600 (Seconds)
 - Manual (selected) / Periodical Key Change
 - Activated From: Sun : 00 : 00
 - Time Interval: 1 (1~168)hour(s)
 - PassPhrase: [Redacted]
 - Confirm PassPhrase: [Redacted]

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D-Link DAP-2660

Home Maintenance Configuration System Logout Help

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